

Laboratorní výsledky pacienta

Datum odběru od - do	11.2.202 6 0 :53	10.2.202 6 0 :20-0 :00	9.2.202 6 05:00	Stav, vývoj	Normály	Jednotky
Metabolismus dusíku						
Urea	5,1		4,2	○ —	3,2 - 8,2	mmol/l
Kreatinin enzymaticky	2		4	● —	5 - 104	μmol/l
eGFR CKD-EPI (kreatinin)	1,910		1,89	○ —	1,5 - 2,5	ml/s/1,73
eGFR (CKD-EPI kreat.) - hodnocení	G1		G1		220 - 547	m ² - μmol/l
Kyselina močová	143		13	● —		mmol/l
Minerály						
Na	4,7		142	○ —	13 - 145	mmol/l
K	104		4,4	○ —	3,5 - 5,1	mmol/l
Cl	8,0		104	○ —	98 - 107	mmol/kg
	3,0					
Osmolalita						
Osmolalita	2,97		282	○ —	280 - 295	μmol/l
Bilirubin						
Bilirubin	1,2			○ —	0 - 25	μmol/l
Bilirubin	4,33			○ —	0 - 4	μkat/l
Bilirubin konjugovaný	0,7, 1,74			○ <		
Enzymy						
ALT	251,9			● >	0,12 -	μkat/l
AST				● —	0,7	μkat/l
GGT				● —	0,22 -	1
ALP				○ —	0 - 1,22	mg/l
FIB4 aELF score				○ —	0,73 -	
FIB-4 skóre				○ —	0 - 1,3	mmol/l
Proteiny						
CRP			2 0,5	●●● —	0 - 10	arb.j.
Glycidy				● —	4,1 - 5,	arb.j.
Glukóza			8,	● —		arb.j.
Hepatitidové markery						
Anti-HAV IgM					4,1 - 5,	arb.j.
Anti-HAV celkové		negativ.			U/l	
HBsAg		pozitiv.			arb.j.	
Anti-HBsAg		negativ.			-	
Anti-HBcAg IgG		2		○	0 - 10	-
Antigen/Anti-HCV		negativ.			arb.j.	
Moč chemicky		později			arb.j.	
pH						
Bílkovina		,0		○	5 -	arb.j.
Krev		1+			arb.j.	
Nitrity-strip		Neg.			arb.j.	
Glukóza		Neg.			arb.j.	
Ketony		Neg.			10 [^] /l	
Urobilinogen		Neg.			10 [^] /l	
Bilirubin		1+			10 [^] /l	
Moč sediment		Neg.			10 [^] /l	
Erytrocyty (M)		0			10 [^] /l	
Leukocyty (M)		50		○	0 - 5 0 -	10 [^] 9/l
Epitelie dlaždicové (M)		20		●	10 0 - 15	10 [^] 12/l
Bakterie		50		●	0 - 40	g/l
Krevní obraz					4 - 10	
Leukocyty přístrojově [WBC]	8,50					
4,17				○	4 - 5,8	
123			12, 0	○ <	135 -	
Erytrocyty přístrojově [RBC]			3,99	○ >		
Hemoglobin [HB]			118	● >	175	

Datum odběru od - do	11.2.202 6 0:53	10.2.202 6 0:20-0:00	9.2.202 6 05:00	Stav, vývoj	Normály	Jednotky
Hematokrit [HCT]	0,383		0,37	● —	0,4 - 0,5	1
Stř.objem erytrocytů [MCV]	91,8		92,0	○ —	82 - 98	fl
Barvivo erytrocytů [MCH]	29,		29,	○ —	28 - 34	pg
Stř.barevná kon.[MCHC]	321		322	○ —	320 - 30	g/l
Erytr.křivka [RDW-CV]	0,132		0,129	○ —	0,1 -	1
Trombocyty přístrojově [PLT] 332			258	○ —	0,152	10 ⁹ /l
Tromb.stř.obj. [MPV]	10,9		11,4	○ <	150 - 400	fl
Trombocytární hematokrit [PCT] 0,3			0,29	● >	7,8 - 11	10 ⁻³
15,9			15,8	○ —	0,12 -	fl
15,9			15,8	○ —	0,35 - 12	10 ⁹ /l
Krev-periferie/příst					18	10 ⁹ /l
Neutrofigly - abs.poč			10,59	○ <<<		10 ⁹ /l
Lymfocyty - abs.poč.			0,83	○ —	2 - 7	10 ⁹ /l
Monocyty - abs.počet			0,9	○ —	0,8 - 4	10 ⁹ /l
Eozinofigly - abs.poč.			0,20	○ —	0,08 - 1,2	10 ⁹ /l
Bazofigly - abs.počet			0,00	○ >	0,001 -	1
Neutrofigly			0,842	● <	0,5	1
Lymfocyt			0,0	● >	0,001 -	1
Monocyt			0,07	○ —	0,2	1
Eozinofigly			0,01	○ —	0,45 - 0,7	1
Bazofigly			0,000	○ —	0,2 - 0,45	-
PCR průkaz DNA/RNA					0,02	-
RNA viru Influenzy A		negativní			0 - 0,05	-
RNA viru Influenzy B		negativní			0 - 0,02	-
COVID-19						
SARS-CoV-2 (PCR) - nosohltan		negativ. *)				arb.j.
Rheumatologie						-
Anti-nukleární IgG /IF/			neděláno *)			arb.j.
Anti-ENA screening			0,5	○	0 - 0,7	kU/l
Anti-dsDNA IgG /IF/			neděláno *)			kU/l
Antiinfekční imunita						kU/l
Anti-Chlamydia pneumoniae			12, *)	●	0 - 11,9	kU/l
Anti-Chlamydia pneumoniae e			24,7	●●●	0 - 11,9	kU/l
Anti-Mycoplasma pneumoniae I gA			18,5	●●●	0 - 11,9	arb.j.
Anti-Mycoplasma pneumoniae I gG			41,4	●●●	0 - 11,9	arb.j.
Anti-Mycoplasma pneumoniae I gM			0,4	○	0 - 0,8	kg/m ³
Nezařazené metody						
Hlen						
Hustota		+				
		1 020				